


Supercross™

PTA Balloon Dilatation Catheter



PRODUCT FEATURES

The flexible tapered tip ensures consistently smooth insertion and atraumatic navigation.

Low crossing profile guarantees excellent crossability to pass through narrow lesions.

Over the wire (OTW) design provides great torque control, kink resistance and pushability to ensure the outstanding deliverability.

Dual-lumen shaft design of the OTW catheter separates the inflation lumen and guide wire lumen to ensure the rapid deflation performance.

Special memory folding technique provides great withdrawal capability, and allows of repeated dilatation.

TECHNICAL PARAMETERS

| | |
|--|---------------------|
| Catheter Design | Rapid Exchange (RX) |
| Balloon Compliance | Semi-compliant |
| Marker | 2 |
| Effective Length | 135cm |
| Minimum Guiding Catheter Compatibility | 5F (0.056") |
| Guide wire Compatibility | 0.014" |

| | |
|--|---------------------|
| Catheter Design | Over the Wire (OTW) |
| Balloon Compliance | Semi-compliant |
| Marker | 2 |
| Effective Length | 80/135cm |
| Minimum Guiding Catheter Compatibility | 5F/6F/7F |
| Guide wire Compatibility | 0.035" |

ORDERING INFORMATION

| Rapid Exchange | | | | | | | |
|-----------------------|---------------------|-----------|-----------|-----------|-----------|-----------|-----------|
| Balloon Diameter (mm) | Balloon Length (mm) | | | | | | |
| | 20 | 30 | 40 | 60 | 80 | 100 | 120 |
| 1.5 | PTAR15020 | / | / | / | / | / | / |
| 2.0 | PTAR20020 | PTAR20030 | PTAR20040 | PTAR20060 | PTAR20080 | PTAR20100 | PTAR20120 |
| 2.5 | PTAR25020 | PTAR25030 | PTAR25040 | PTAR25060 | PTAR25080 | PTAR25100 | PTAR25120 |
| 3.0 | PTAR30020 | PTAR30030 | PTAR30040 | PTAR30060 | PTAR30080 | PTAR30100 | PTAR30120 |
| 3.5 | PTAR35020 | PTAR35030 | PTAR35040 | PTAR35060 | PTAR35080 | PTAR35100 | PTAR35120 |
| 4.0 | PTAR40020 | PTAR40030 | PTAR40040 | PTAR40060 | PTAR40080 | PTAR40100 | PTAR40120 |
| 4.5 | PTAR45020 | PTAR45030 | PTAR45040 | PTAR45060 | PTAR45080 | PTAR45100 | PTAR45120 |
| 5.0 | PTAR50020 | PTAR50030 | PTAR50040 | PTAR50060 | PTAR50080 | PTAR50100 | PTAR50120 |
| 5.5 | PTAR55020 | PTAR55030 | PTAR55040 | PTAR55060 | PTAR55080 | PTAR55100 | PTAR55120 |
| 6.0 | PTAR60020 | PTAR60030 | PTAR60040 | PTAR60060 | PTAR60080 | PTAR60100 | PTAR60120 |
| 7.0 | PTAR70020 | PTAR70030 | PTAR70040 | PTAR70060 | PTAR70080 | PTAR70100 | PTAR70120 |

| Over the Wire (Effective length is 80cm) | | | | | | | |
|--|---------------------|-----------|-----------|-----------|-----------|-----------|-----------|
| Balloon Diameter (mm) | Balloon Length (mm) | | | | | | |
| | 20 | 30 | 40 | 60 | 80 | 100 | 120 |
| 3.0 | PTAO30020 | PTAO30030 | PTAO30040 | PTAO30060 | PTAO30080 | PTAO30100 | PTAO30120 |
| 3.5 | PTAO35020 | PTAO35030 | PTAO35040 | PTAO35060 | PTAO35080 | PTAO35100 | PTAO35120 |
| 4.0 | PTAO40020 | PTAO40030 | PTAO40040 | PTAO40060 | PTAO40080 | PTAO40100 | PTAO40120 |
| 4.5 | PTAO45020 | PTAO45030 | PTAO45040 | PTAO45060 | PTAO45080 | PTAO45100 | PTAO45120 |
| 5.0 | PTAO50020 | PTAO50030 | PTAO50040 | PTAO50060 | PTAO50080 | PTAO50100 | PTAO50120 |
| 5.5 | PTAO55020 | PTAO55030 | PTAO55040 | PTAO55060 | PTAO55080 | PTAO55100 | PTAO55120 |
| 6.0 | PTAO60020 | PTAO60030 | PTAO60040 | PTAO60060 | PTAO60080 | PTAO60100 | PTAO60120 |
| 7.0 | PTAO70020 | PTAO70030 | PTAO70040 | PTAO70060 | PTAO70080 | PTAO70100 | PTAO70120 |
| 8.0 | PTAO80020 | PTAO80030 | PTAO80040 | PTAO80060 | PTAO80080 | PTAO80100 | PTAO80120 |
| 9.0 | PTAO90020 | PTAO90030 | PTAO90040 | PTAO90060 | PTAO90080 | PTAO90100 | PTAO90120 |
| 10.0 | PTAO10020 | PTAO10030 | PTAO10040 | PTAO10060 | PTAO10080 | PTAO10100 | PTAO10120 |
| 12.0 | PTAO12020 | / | / | / | / | / | / |

| Over the Wire (Effective length is 135cm) | | | | | | | |
|---|---------------------|-------------|-------------|-------------|-------------|-------------|-------------|
| Balloon Diameter (mm) | Balloon Length (mm) | | | | | | |
| | 20 | 30 | 40 | 60 | 80 | 100 | 120 |
| 3.0 | PTAO30020 L | PTAO30030 L | PTAO30040 L | PTAO30060 L | PTAO30080 L | PTAO30100 L | PTAO30120 L |
| 3.5 | PTAO35020 L | PTAO35030 L | PTAO35040 L | PTAO35060 L | PTAO35080 L | PTAO35100 L | PTAO35120 L |
| 4.0 | PTAO40020 L | PTAO40030 L | PTAO40040 L | PTAO40060 L | PTAO40080 L | PTAO40100 L | PTAO40120 L |
| 4.5 | PTAO45020 L | PTAO45030 L | PTAO45040 L | PTAO45060 L | PTAO45080 L | PTAO45100 L | PTAO45120 L |
| 5.0 | PTAO50020 L | PTAO50030 L | PTAO50040 L | PTAO50060 L | PTAO50080 L | PTAO50100 L | PTAO50120 L |
| 5.5 | PTAO55020 L | PTAO55030 L | PTAO55040 L | PTAO55060 L | PTAO55080 L | PTAO55100 L | PTAO55120 L |
| 6.0 | PTAO60020 L | PTAO60030 L | PTAO60040 L | PTAO60060 L | PTAO60080 L | PTAO60100 L | PTAO60120 L |
| 7.0 | PTAO70020 L | PTAO70030 L | PTAO70040 L | PTAO70060 L | PTAO70080 L | PTAO70100 L | PTAO70120 L |
| 8.0 | PTAO80020 L | PTAO80030 L | PTAO80040 L | PTAO80060 L | PTAO80080 L | PTAO80100 L | PTAO80120 L |
| 9.0 | PTAO90020 L | PTAO90030 L | PTAO90040 L | PTAO90060 L | PTAO90080 L | PTAO90100 L | PTAO90120 L |
| 10.0 | PTAO10020 L | PTAO10030 L | PTAO10040 L | PTAO10060 L | PTAO10080 L | PTAO10100 L | PTAO10120 L |
| 12.0 | PTAO12020 L | / | / | / | / | / | / |

